

Applications

- Disposal of contaminated sharps

Features

- Designed and labelled in accordance to the UN3291 standard
- Tested according to ISO 23907:2019
- Distinctive colouring yellow polypropylene container and clear labelling
- Accommodates a wide variety of sizes of sharps.
- Stable flat surface
- Resistant to solvents and perforation
- Easy to handle and operate rotating open/close
- Leak-proof upon locking the lid and puncture-resistant






Aspects

- Box: Yellow/Red
- Lid: Yellow/Transparent
- Button: Yellow/Red
- Ring: Blue

Packaging Specifications

Product Code	Product Description	Containers/ Box	Covers/ Box	Cardboard Box Dimensions (mm)	Average Weight/Box (Kg)
SH-100	Sharps Containers 10L	20	20	454 x 454 x 593	10.75
SH-050 SH-050R	Sharps Containers 5L	30	30	689 x 499 x 398	12.14
SH-030	Sharps Containers 3L	30	30	513 x 512 x 400	11.07
SH-018	Sharps Containers 1.8L	45	45	513 x 512 x 400	11.34
SH-010	Sharps Containers 1L	100	100	435 x 435 x 49.5	12.77

Technical Specifications

Product	Average Weight (g)	Length	Width	Height
 <p>Sharp container 10L (SH-100)</p>	481	230.0 (±0.5)	225.0 (±0.5)	305.0 (±0.5)
 <p>Sharp container 5L (SH-050)</p>	366	230.0 (±0.5)	225.0 (±0.5)	190.0 (±0.5)
 <p>Sharp container 3L (SH-030)</p>	281	174.0 (±0.5)	167.0 (±0.5)	233.0 (±0.5)
 <p>Sharp container 1.8L (SH-018)</p>	233	174.0 (±0.5)	167.0 (±0.5)	163.0 (±0.5)
 <p>Sharp container 1L (SH-010)</p>	114.3	90.0 (±0.5)	90.0 (±0.5)	195.0 (±0.5)