



Applications

- For Collection, transportation, and storage of a wide variety of specimens.

Features

- Disposable single use containers
- Polypropylene or polystyrene container with screw cap
- Inert and chemically resistant to commonly used laboratory reagents at room temperature
- Label for patient's information Legibility reducing errors in the preanalytical phase
- Provide an integrity seal preventing leakage problems
- Sterile containers are individually wrapped with BOPP film

Product and Packaging Specifications

Product Code	Material Type	Sterility	Spoon	Label	Cap color	Pcs/ Bag	Bags /Box	Pcs/ Box	Box Dimensions L x w x H (mm)	Average Weight/Box (Kg)		
SCP-600	Polypropylene	Aseptic	No	No	Red	500	1	500	414 x 414 x 346 (± 10)	6.7 (± 0.5)		
SCP-610					White							
SCP-640				Yes	Red							
SCP-650					White							
SCP-601		Sterile EO		Yes	No	Red	1				500	
SCP-611						White						
SCP-641					Yes	Red						
SCP-651					Yes	White						
SCP-620		Aseptic	Yes		No	Red	500				1	
SCP-630						White						
SCP-670						Yes						White
SCP-621						Sterile EO						No
SCP-631		White										
SCP-661		Yes		Red								
SCP-661		Yes		Red								
SCS-600	Polystyrene	Aseptic		No	No	Red	500	1	500	414 x 414 x 346 (± 10)	6.7 (± 0.5)	
SCS-610			White									
SCS-640			Yes		Red							
SCS-650					White							
SCS-601		Sterile EO	Yes		No	Red	1	500				
SCS-611						White						
SCS-641					Yes	Red						
SCS-651					Yes	White						
SCS-620		Aseptic		Yes	No	Red	500	1				
SCS-630						White						
SCS-660						Yes						Red
SCS-661						Sterile EO						Yes
SCS-671		White										