



### Applications

- For Collection, transportation, and storage of a wide variety of specimens.

### Features

- Disposable containers
- Molded in clear polypropylene
- Inert and chemically resistant to commonly used laboratory reagents at room temperature
- Legibility of the patient's information on the label reducing errors in the preanalytical phase
- Provide an integrity seal preventing leakage problems

### Product and Packaging Specifications

Product Code	Size	Sterility	Spoon	Label	Cap color	Pcs /Bag	Bags /Box	Pcs /Box	Box Dimensions : L x w x H (mm)	Box Average Weight (Kg)				
SSP-630	60 ml	Aseptic	White spoon	No	White cap	500	1	500	414 x 414 x 346 (± 10)	6.7 (±0.5)				
SSP-631		Sterile				1	500							
SSP-670		Aseptic		500		1								
SSP-671		Sterile		1		500								
SSP-130	120 ml	Aseptic		White spoon		No	White cap	500	1	500	512 x 500 x 280 (±10)	10 (±0.5)		
SSP-131		Sterile						1	250	250	599 x 299 x 408 (±10)	5.2 (±0.5)		
SSP-530	150 ml	Aseptic						No	White cap	500	1	500	564 x 424 x 295 (±10)	12 (±0.5)
SSP-531		Sterile								1	250	250	599 x 299 x 495 (±10)	6.4 (±0.5)

### Technical Specifications