



Applications

- For Collection, transportation, and storage of a wide variety of specimens.

Features

- Disposable containers
- Molded in clear polypropylene container with screw cap
- Inert and chemically resistant to commonly used laboratory reagents at room temperature
- Label for patient's information Legibility reducing errors in the preanalytical phase
- Provide an integrity seal preventing leakage problems
- Sterile containers are individually wrapped with BOPP film

Product and Packaging Specifications

Product Code	Sterility	Spoon	Label	Cap color	Pcs/Bag	Bags/Box	Pcs/Box	Box Dimensions L x w x H (mm)	Average Weight/Box (Kg)
SCP-100	Aseptic	Without spoon	No	Red	Bulk 500 containers / 500 caps	2	500	512 x 500 x 280 (±10)	7.35 (±0.5)
SCP-110				White					
SCP-101	Sterile			Red	Individually wrapped	250	250	599 x 299 x 408 (±10)	3.84 (±0.5)
SCP-111				White					
SCP-120	Aseptic	With spoon		Red	Bulk 500 containers / 500 caps	2	500	512 x 500 x 280 (±10)	7.35 (±0.5)
SCP-130				White					
SCP-121	Sterile			Red	Individually wrapped	250	250	599 x 299 x 408 (±10)	3.84 (±0.5)
SCP-131				White					
SCP-141		Yes	Red						