



Applications

- For Collection, transportation, and storage of a wide variety of specimens.

Features

- Disposable single use containers
- Molded in clear polypropylene containers with screw cap
- Frosted Label
- Inert and chemically resistant to commonly used laboratory reagents at room temperature
- Provide an integrity seal preventing leakage problems
- Sterile containers are individually wrapped with BOPP film

Product and Packaging Specifications

Product Code	Sterility	Spoon Availability	Spoon color	Label	Cap color	Pcs /Bag	Bags /Box	Pcs /Box	Box Dimensions L x W x H (mm)	Average Weight/Box (Kg)		
SCP-500	Aseptic	No	NA	No	Red	500	1	500	564 x 424 x 295 (±10)	12 (±0.5)		
SCP-510					White							
SCP-501	Sterile				Red	1	250	250	599 x 299 x 495 (±10)	6.4 (±0.5)		
SCP-511					White							
SCP-520	Aseptic	Yes	Red		No	Red	500	1	500	564 x 424 x 295 (±10)	12 (±0.5)	
SCP-530			White			White						
SCP-521	Sterile		Red			No	Red	1	250	250	599 x 299 x 495 (±10)	6.4 (±0.5)
SCP-531			White				White					
SCP-541			NA	Yes	Red							