


MICROPIPETETTE TIPS**Applications**

- Used in the extraction of small and accurate amounts of liquid.

Features

- Precision molded tips allowing complete retrieval of precious samples
- Consistent volume accuracy
- Packaged either in bulk in economical PE bags or in racks
- Match a variety of pipettes
- Flash free orifice avoiding nucleic acid and protein denaturation
- Yellow and blue standard colors are for easy identification. Other colors can be manufactured upon request
- Can be supplied with Pipette Filters made of HDPE & PP to help prevent liquid bypass and airborne contaminants from being drawn into the tips.

Product and Packaging Specifications
TIPS:

Product Code	Color	Volume (µL)	Filter	Sterility	Packaging	Pcs /Bag	Bags /Box	Pcs /Box	Box Dimensions L x w x H (mm)	Average Weight/ Box (Kg)
MT-101	Blue	1000	No	No	Bulk in PE bag	1000	16	16000	604 x 404 x 398	11.52
MT-111			Yes	Yes						13.13
MT-201			No	No	In rack	96	30	2880	340 x 335 x 440	5.47
MT-211			Yes	Yes						5.76
MT-102	Yellow	5-200	No	No	Bulk in PE bag	1000	30	30 000	584 x 334 x 318	11.66
MT-112		200	Yes	Yes						12.57
MT-122		100	Yes	Yes						12.52
MT-202		5-200	No	No	In rack	96	50	4800	454 x 254 x 358	5.58
MT-212		200	Yes	Yes						5.73
MT-222		100	Yes	Yes						5.68

RACKS:

Product code	Color	Item	Pcs/box	Box Dimensions L x w x H (mm)	Average Weight/ Box (Kg)
MT-301	Blue	Rack For Blue Tips	54	340 x 335 x 440	3.59
MT-302	Yellow	Rack For Yellow Tips	96	454 x 254 x 358	3.91