



Applications

- Plastic single use loops for inoculation of colonies to growth media or for staining slides

Features

- Acid Resistant
- Choice of hard/rigid or soft/flexible loops to cover all microbiologist applications
- Unique design providing a consistent volume transfer
- Uniform and smooth streaking without agar cutting
- Designed for a double use: one side for transporting colonies and the other for petri dishes streaking
- Loops' shaft is designed with hexagonal shape for better grip and easy control
- Packed in sterile peel pouch
- Supplied in different packaged quantity per pouch (1, 5, 10 or 20 pcs/pouch) or pack.
- Calibration certificate can be provided for accuracy and precision using Evans Blue Dye Method

Product and Packaging Specifications

| Design | Color | Sterility | Product code | Pcs /Pouch | Pouches /Packet | Packets /Box | Pcs /Box | Box Dimensions L x H x W (mm) | Avrg Box Weight (Kg) |
|-------------|--|-----------|--------------|------------|-----------------|--------------|----------|--------------------------------|----------------------|
| Soft | Yellow/ customer requirement | Sterile | IL-011 | 1 | 250 | 5 | 1250 | 414 x 314 x 288 (\pm 10) | 5 |
| | | | IL-0111 | 1 | 500 | 5 | 2500 | | 5.5 |
| | | | IL-012 | 5 | 200 | 5 | 5000 | | 8 |
| | | | IL-013 | 10 | 100 | 5 | | | 7 |
| | | | IL-014 | 20 | 50 | 5 | | | 6.6 |
| Hard | Golden Yellow/ customer requirement | Sterile | IL-015 | 1 | 250 | 5 | 1250 | | 5 |
| | | | IL-016 | 5 | 200 | 5 | 5000 | | 8 |
| | | | IL-017 | 10 | 100 | 5 | | | 7 |
| | | | IL-018 | 20 | 50 | 5 | | | 6.6 |