

MEDICINE CUP**Applications**

- Used for oral medicines administration

Features

- Manufactured from polypropylene (PP)
- Graduated, 30ml capacity
- Pliable, unbreakable
- Caps available separately
- Transparent is the standard color, other colors available on order

Product and Packaging Specifications

Product Code	Product Description	Pcs/Bag	Bags/Box	Pcs/Box	Carton Dimensions(mm)	Average Weight/Box(kg)
MC-100	Medicine cup without lid	50	100	5000	405 x 310 x 380 (± 10)	12.5
MC-045	Medicine cup without lid	50	90	4500	405 x 335 x 370 (± 10)	11.5
MC-101	Lid for medicine cup	10	100	1000	303 x 270 x 85 (± 10)	1.27
MC-102	Lid for medicine cup	25	100	2500	330 x 250 x 260 (± 10)	3.11