



Applications

- Used for the constriction and control of blood circulation

Features

- Ergonomic
- Easy to use
- Blue tape with white buckle

Material Type(s)

- Tape: Cotton & Terylene Silk, Latex Free Rubber
- Buckle: ABS

Product and Packaging Specifications

Product Code	Weight(g)	Pcs/Bag	Bags/Box	Pcs/Box	Carton Dimensions (mm)	Average Weight/Box
VBCTQ-100	30.5	1	300	300	588 x 400 x 200	11.2 Kg