

UMBILICAL CORD CLAMP



Applications

- Used in clamping the umbilical of newly born babies
- It can be used as a disposable first aid item and easily wiped down with antiseptic.

Features

- White is standard color, other colors available on order
- Acetal copolymer material.
- Single use.
- Tamper-proof security lock to prevent accidental opening after clamping
- Individually wrapped
- Disposable
- Autoclavable
- Graduated
- Flattened profile ensuring comfortable carrying
- Different opaque colors are available
- Water proof
- Heat resistant

Product and Packaging Specifications

Product Code	Sterility	Material Type	Pcs/pack	packs/Box	Pcs/Box	Carton Dimensions (mm)	Average Weight/Box (Kg)
CC-100	ETO	Acetal	75 Individually wrapped pcs.	40	3000	575 x 510 x 425 (±10)	13.5