



Applications

- Tubes for blood collection.

Features

- Transparent Polypropylene tubes
- Push Polyethylene cap for a leak free closure
- For 2.5 ml or 5ml volume of blood.
- Color coded label and cap for ease of recognition
- Round bottom
- Supplied plain or with an anticoagulant
- CE marked in accordance with (EU) 2017/746 2017
- Supplied with rack for user convenience

Product and Packaging Specifications

Product code	Volume on label (ml)	Closure color	Anticoagulant	Application	Pcs/ Rack	Pcs/ Box	Box Dimensions L x w x H (mm)	Box Average Weight (Kg)				
BCP-900	No label	Grey	PLAIN	NA	50	1000	464 x 179 x 368 (±10)	4.25 (±0.5)				
BCP-210	2.5 / 5											
BCP-211	2.5	Yellow	FLUORIDE OXALATE	Glucose determinations on Plasma.								
BCP-511	5											
BCP-212	2.5	Green or Violet	K3EDTA	For whole blood Hematology determinations								
BCP-512	5											
BCP-213	2.5	Dark Blue	LITHIUM HEPARIN	for Plasma determinations					50	1000	464 x 179 x 368 (±10)	4.25 (±0.5)
BCP-513	5											
BCP-214	2.5	Pink	TRISODIUM CITRATE	Routine coagulation studies								
BCP-514	5.0											